

MINIMISED TRAUMA & ENHANCED PATIENT COMFORT


Designed for Initial Endoscopic, Radiologic
or Laparoscopic G, J & GJ Tube Placement (12-20 Fr)



Introducer Kits provide reduced surgical steps and fewer incisions^{1,2,3}

SAF-T-PEXY* NEEDLES:


- **Simple, Intuitive Tools**
 - Confidently secure the stomach
 - Audible confirmation of suture lock
 - Suture lock can bidirectionally be adjusted during application
- **Promotes Patient Compliance & Comfort**
 - Absorbable sutures
 - Promotes patient hygiene with easy-to-clean bolsters
 - Smooth, round retention device for patient comfort
- **Controlled Delivery**
 - One-handed, on-demand T-bar deployment
 - Pre-loaded needles: protects clinician from inadvertent needle sticks



Lowered risk for surgical & post-op complications⁴

SERIAL DILATOR:

- Allows Seldinger Technique:
Reduced need for repeat passing
over the guide wire
- Available in 5 different French sizes:
Fr 16, 18, 20, 22, 24
Fewer tools to handle
Reach desired dilation with a
single telescoping device
- Integrated Peel Away Sheath:
Facilitates feeding tube placement
- Made to fit a wide range of tubes:
MIC* & MIC-KEY* G, J or GJ tubes
Fr 12 up to Fr 20



AVANOS* MIC-KEY* OVER-THE-WIRE STOMA MEASURING DEVICE

NAME	PACKAGING	PRODUCT CODE
MIC-KEY* Stoma Measuring Device	10 Each/Kit	98460




- Ensures stoma tract length accuracy
- Dual lumen construction allows Seldinger technique
- Proper fit results in long-lasting patient comfort



Delivers all the essential tools for an efficient, secure and less invasive placement

KIT OFFERS:

- Lowered risk for surgical and post-op complications than other procedures (e.g. longer hospital stay, punctured colon, tissue damage, etc.)⁴
- Tools specifically designed for visual and surgical accuracy: SAF-T-PEXY* T-Fasteners, Gastrointestinal Anchor Set, Telescopic Dilator with integrated Peel Away Sheath
- Sterilised introducer needle, hemostat, J-guide wire/seeking catheter, scalpel and more
- Flexibility as replacement or primary procedure of a wide range of tubes: MIC* & MIC-KEY* G, J or GJ Tube, Fr 12 up to Fr 20^{5,6}
- Easy access to the essential tools in one convenient kit



AVANOS* INTRODUCER KITS FOR GASTROSTOMY FEEDING TUBES

NAME	PACKAGING	PRODUCT CODE
Introducer Kit for Gastrostomy Feeding Tube, 16 Fr Dilator	1 kit/case	98430
Introducer Kit for Gastrostomy Feeding Tube, 18 Fr Dilator	1 kit/case	98431
Introducer Kit for Gastrostomy Feeding Tube, 20 Fr Dilator	1 kit/case	98432
Introducer Kit for Gastrostomy Feeding Tube, 22 Fr Dilator	1 kit/case	98433
Introducer Kit for Gastrostomy Feeding Tube, 24 Fr Dilator	1 kit/case	98434

NAME	PACKAGING	PRODUCT CODE
Laparoscopic Introducer Kit for Gastrostomy Feeding Tube, 16 Fr Dilator	1 kit/case	10250
Laparoscopic Introducer Kit for Gastrostomy Feeding Tube, 18 Fr Dilator	1 kit/case	10230
Laparoscopic Introducer Kit for Gastrostomy Feeding Tube, 20 Fr Dilator	1 kit/case	10240
Laparoscopic Introducer Kit for Gastrostomy Feeding Tube, 22 Fr Dilator	1 kit/case	10237
Laparoscopic Introducer Kit for Gastrostomy Feeding Tube, 24 Fr Dilator	1 kit/case	10226

AVANOS* INTRODUCER KITS FOR JEJUNAL/GASTRIC JEJUNAL FEEDING TUBES

NAME	PACKAGING	PRODUCT CODE
Introducer Kit for Jejunal/Gastric Jejunal Feeding Tube, 18 Fr Dilator	1 kit/case	98436A
Introducer Kit for Jejunal/Gastric Jejunal Feeding Tube, 20 Fr Dilator	1 kit/case	98437
Introducer Kit for Jejunal/Gastric Jejunal Feeding Tube, 22 Fr Dilator	1 kit/case	98438

NAME	PACKAGING	PRODUCT CODE
Laparoscopic Introducer Kit for Jejunal/Gastric Jejunal Feeding Tube, 18 Fr Dilator	1 kit/case	10269A
Laparoscopic Introducer Kit for Jejunal/Gastric Jejunal Feeding Tube, 20 Fr Dilator	1 kit/case	10249
Laparoscopic Introducer Kit for Jejunal/Gastric Jejunal Feeding Tube, 22 Fr Dilator	1 kit/case	10229

AVANOS* SAF-T-PEXY* T-FASTENERS

NAME	PACKAGING	PRODUCT CODE
SAF-T-PEXY* T-Fasteners, Gastrointestinal Anchor Set	1 each (3 needles + 1 curved hemostat)	98701

AVANOS* SERIAL DILATORS

NAME	PACKAGING	PRODUCT CODE
Serial Dilator with Built-in Peel-Away Sheath, 16 Fr	1 each	98702
Serial Dilator with Built-in Peel-Away Sheath, 18 Fr	1 each	98703
Serial Dilator with Built-in Peel-Away Sheath, 20 Fr	1 each	98704
Serial Dilator with Built-in Peel-Away Sheath, 22 Fr	1 each	98705
Serial Dilator with Built-in Peel-Away Sheath, 24 Fr	1 each	98706



1. Michaud L., et al.: One-step Percutaneous Gastrojejunostomy in Early Infancy. JPGN, 2012, vol. 54, pp. 820-821.
2. Gothberg G, et. Al.: One-Step Insertion of Low-Profile Gastrostomy in Pediatric Patients vs Pull PEG: Retrospective Analysis of outcomes. Journal of Parenteral and Enteral Nutrition, JPEN, February, 2015.
3. When compared to COOK® Medical enteral feeding tube placement kits and devices
4. Antonoff MB, Hess DJ, Saltzman DA, Acton RD. Modified approach to laparoscopic gastrostomy tube placement minimizes complications. Pediatr Surg Int 2009. Apr;25(4):349-353.10.1007/s00383-009-2340-z
5. Michaud L., et al.: One-step Percutaneous Gastrojejunostomy in Early Infancy. JPGN, 2012, vol. 54, pp. 820-821
6. Gothberg G, et. Al.: One-Step Insertion of Low-Profile Gastrostomy in Pediatric Patients vs Pull PEG: Retrospective Analysis of outcomes. Journal of Parenteral and Enteral Nutrition, JPEN, February, 2015.

**From the UK, please call Avanos
at 0800 917 65 85 - fax 0800 169 02 35
or send an email to
customerservice.uk.ie@avanos.com**

**From IRL, please call Uniphar Hospital Services
at +353 1 468 84 56 - fax +353 1 466 07 82
or send an email to unipharhospitalervices@uniphar.ie**

**From all other countries, please call Avanos
at +32 2 700 68 51 - fax +32 2 711 26 91
or send an email to customerservice.export@avanos.com**

AVANOS

**www.mic-key.com
www.avanos.co.uk**

* Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2019 AVNS. All rights reserved.

All other trademarks are the property of their respective owners. Biosyn™ is a Registered Trademark of United States Surgical Corporation. HC218-00-UK